



# Tia Dent<sup>®</sup>



## Tia Dent

Special dental arch stone

Type 4 according to DIN EN ISO 6873

Discover the new dentale arch stone with fascinating properties.  
For extremely fine processing - with maximum hardness as a result.

Mixing ratio for 100 g powder:	20 - 25 ml distilled water
Working time:	6 - 7 minutes
End of setting:	approx. 35-40 minutes
Hardness after 24 h:	approx. 330 MPa (N/mm <sup>2</sup> )
Compression strength:	approx. 70 MPa (N/mm <sup>2</sup> )
Expansion acc. to DIN EN (2 h):	< 0.08 %



For More Information Please Contact us at 0098 912 4000 739 ( phone and Whatsapp )

 [Safavigypsum.ir/tiadent](http://Safavigypsum.ir/tiadent)

**Tia Dent**